

Upgrade your hearing

Choose the Interton Crisp that suit your needs

TYPICAL LISTENING SITUATIONS	CRISP ²	CRISP ³
Speech in quiet surroundings	● ● ● ● ● ● ● ●	● ● ● ● ● ● ● ●
Situations without speech	● ● ● ● ● ● ● ●	● ● ● ● ● ● ● ●
Listening to music	● ● ● ● ● ● ● ●	● ● ● ● ● ● ● ●
Closeness to person or object	● ● ● ● ● ● ● ○	● ● ● ● ● ● ● ●
Wireless accessories	● ● ○ ○ ○ ○ ○ ○	● ● ● ● ● ● ● ●
Speech in background noise	● ● ● ● ● ● ○ ○	● ● ● ● ● ● ● ○
Office situations	● ● ● ● ○ ○ ○ ○	● ● ● ● ● ● ○ ○
Changing conversation focus	● ● ● ● ○ ○ ○ ○	● ● ● ● ● ● ○ ○
Meeting situations	● ● ● ● ○ ○ ○ ○	● ● ● ● ● ● ○ ○
Personalization options	● ● ● ● ○ ○ ○ ○	● ● ● ● ● ● ○ ○

Interton Crisp improves everyday communication with better speech understanding in noisy situations and no feedback issues. With Crisp 3 you also get access to all the wireless accessories and the small Receiver-in-the-ear.

Interton - Essential Hearing

Founded in Germany in 1962, Interton today is part of one of the world's largest hearing instrument organizations, the GN Group. We use their proven core technology in a way that falls in line with our own way of doing things. We believe that most people need a straight-forward, functional and easy-to-use hearing solution at a fair price. So that's what we give them.

Interton A/S | Lautrupbjerg 7 | DK-2750 Ballerup
 Tel: +45 4575 1111 | www.interton.com

M101416-GB - INTERTON AND INTERTON CRISP ARE REGISTERED TRADEMARKS - RAFFINADERIET.DK

Everyone deserves an upgrade!



Every year we all lose a small part of our hearing, and the rate of loss increases with age. If conversations are getting hard to follow, maybe it's time to upgrade your hearing?

It could be time to invest in a hearing aid that's discreet and easy to use. One that makes speech clearer and conversation easier, but is almost unnoticeable. One with wireless connectivity to help you enjoy the many electronic devices in your life without disturbing others.

Interton Crisp gives you all these benefits you normally find in much more expensive hearing aids, but at a price you can afford! Better hearing, plus satisfying value for money. That's the kind of upgrade we all appreciate.



Watch the Interton Crisp movie

All you need is to download the free Scanlife® application and simply scan the code.

You can also watch the movie at www.interton.com

Great sound even in noisy surroundings

Interton's innovative technologies make hearing in everyday situations so much easier. Crisp brings back lost high frequencies and filters out background noise. It restores speech clarity, while maintaining a relaxed, natural sound, so you can follow a conversation even in noisy places like restaurants or airports.

SoundScape technology brings back the natural sound experience. It mimics the organic sound processing of the human ear to produce sounds that are as rich, clear and natural as they can be. That means you hear everything that's going on around you, and can locate where sounds are coming from. Our Feedback Management System also eliminates disturbing whistling sounds when you put things close to your ears or hug your loved ones.



Better hearing Almost invisible

Hearing loss is a natural part of aging, but it doesn't mean you have to stop doing the things you enjoy. Using a tiny, discreet Crisp model, you can get all the benefits of quality sound and easy social interaction.



Crisp lets you choose the type of hearing aid you are most comfortable with. Whether you prefer a custom-made in the ear model or one that sits neatly and discreetly behind the ear, a wide range of color options is available.

Choose a solution which matches your skin or hair color perfectly making the hearing aid practically invisible.



Colors:



Beige

*Medium
Blonde*



*Dark
Brown*

*Marble
Grey*



Black

Wet weather? No problem!

Crisp is tested under all kinds of conditions to simulate the rough and tumble of real life, and it comes through every time. What makes Crisp so hardwearing?

Key parts of the hearing aid are protected by a shell reinforced by glass fibers which makes it extremely tough. What's more, all our hearing aids are treated with a thin layer of Nano coating. This shields every part of the hearing aid - including the electronics on both the inside and the outside - so when humidity or sweat come into contact with it, they form beads and just roll right off.



Great wireless opportunities

Crisp's wireless connectivity gives you total freedom and control, with a range of easy-to-use accessories. The handy Remote Control lets you manage your hearing aids easily by adjusting the volume, changing programs or switching between your TV, radio or computer.

The TV Streamer lets you enjoy TV, stereo or audio from your computer without disturbing others. Enjoy high quality crisp stereo sound directly in your ears, without wearing uncomfortable adapters around your neck. You will still be able to participate in conversations around you.

A tiny Phone Clip connects your Crisp to a Bluetooth®-enabled mobile phone, allowing you to answer or end calls as well as change the volume at the touch of a button.



Interton
Remote Control

Interton
TV Streamer

Interton
Phone Clip